

**mindray**

## **CL-2000i**

Chemiluminescence Immunoassay System



# Mindray Solution for Immunoassay Tests

01 >>

## High Throughput

Up to 240 tests/hour

04 >>

## Flexible Scalability

Flexible connection with BS-800 module and BS-2000 module into integrated system



02 >>

## Large Capacity

- 300 samples in one batch
- Support continuous loading
- 36 onboard capacity for reagent cartridges
- 1134 cuvettes loaded in one batch

03 >>

## Easy Operation

- Intuitive interface for software
- Real-time indication of cuvettes, waste and wash buffer status
- Step-by-step maintenance guide


# Test menu

|   |   |  |   |  |   |
|---|---|--|---|--|---|
| Thyroid   |                                | Fertility  |    | Bone Metabolic   |    |
|    | <b>Infectious Disease</b>   |   | <b>Tumor marker</b>   |   | <b>Cardiac</b>  |
| FT3<br>FT4<br>T3<br>T4<br>TSH<br>Anti-Tg<br>Anti-TPO<br>Tg<br>rT3<br>TRAb           | HIV Combo<br>HBsAg<br>Anti-HBs<br>HBeAg<br>Anti-HBe<br>Anti-HBc<br>Anti-HCV**<br>Anti-TP (Syphilis)<br>HAV IgM* | Total $\beta$ -HCG<br>FSH<br>LH<br>Prolactin<br>Estradiol<br>Estriol<br>TESTO<br>PROG<br>AMH<br>Free testosterone*<br>17-OH PROG*<br>SHBG* | CEA<br>AFP<br>CA125<br>CA15-3<br>CA19-9<br>Free PSA<br>Total PSA<br>NSE<br>CYFRA 21-1<br>CA72-4<br>PG I<br>PG II<br>SCCA<br>HE4<br>ProGRP<br>CA50<br>CA242<br>PIVKA-II*<br>AFP-L3%*<br>Hp IgG*<br>G-17* | Intact PTH<br>Calcitonin<br>Vitamin D total  | CK-MB<br>Myoglobin<br>Troponin I<br>BNP<br>NT-proBNP*                                 |
| Diabetes  |                              | Liver Fibrosis   |    | ToRCH  |  |
|  | <b>Adrenal Function</b>   |   | <b>COVID-19</b>   |   | <b>Growth Hormone</b>   |
| DM<br>Insulin<br>C-peptide  | DHEA-S<br>Cortisol<br>ACTH  | Laminin<br>Hyaluronic Acid<br>PIIINP<br>Collagen IV  | SARS-COV-2 IgG<br>SARS-COV-2 IgM<br>SARS-COV-2S-RBD IgG<br>SARS-COV-2 Neutralizing Antibody   | Toxo IgG<br>Toxo IgM<br>Rubella IgG<br>Rubella IgM<br>CMV IgG<br>CMV IgM<br>HSV-1/2 IgG*<br>HSV-1/2 IgM*<br>HSV-1 IgG*<br>HSV-2 IgG* | GH*<br>IGF-1*   |
| Anemia  |                              | Inflammation   |    |  |   |
|  | <b>Hypertension</b>   |   | <b>DS' Screening</b>  |  |   |
| Ferritin<br>Vitamin B12<br>Folate<br>RBC Folate                                     | Renin<br>Aldosterone  | Procalcitonin<br>IL-6<br>sCD14-ST  | PAPP-A*<br>Free $\beta$ -HCG*   |  |   |

\*: Under development  
\*\*: Non-CE

# CL-2000i

## Chemiluminescence Immunoassay System



### Manipulator 2

Moving cuvettes for mixing, magnetic separation and incubation

### Reaction Area

100 positions for incubation  
50 positions for reaction and reading  
Direct solid heating

### Reagent Probe

20 $\mu$ l-135 $\mu$ l, with increment of 1 $\mu$ l  
Bubble detection, liquid level detection  
Collision protection and auto recovery  
Enhanced exterior and interior washing

### Reagent Carousel

36 onboard capacity for reagent cartridges  
2 $^{\circ}$ C~8 $^{\circ}$ C constant cooling compartment  
Continuous reagent loading

### Normal Lane

### Return Lane

### Passing Lane

Higher priority for STAT, calibrator, control and rerun racks  
For routine sample racks transferred to other analytical unit(s)

### Sample Probe

10 $\mu$ l-200 $\mu$ l, with increment of 1 $\mu$ l  
Clot detection, liquid level detection and tracking  
Collision protection and auto recovery  
Enhanced exterior and interior washing

### Cuvettes Area

9 boxes of cuvettes can be loaded in one batch  
126 cuvettes per box  
Automatic unloading for empty boxes

# System Layout

## Mixing Unit

High speed vortexer mixing for reaction solution

## Magnetic Separation Unit

4-phase washing

## Sample Delivery Module (SDM)

Up to 300 samples to be loaded  
Continuous sample loading and offloading

## STAT and RUN buttons

Flexible STAT rack loading  
One touch to initiate analysis

## Substrate Bottles

2 bottles for equilibrium  
2 bottles ready for aspirating  
Bottom piercing to eliminate contamination

## Manipulator 1

Moving cuvettes for dispensing  
Intelligent cuvette detection




# Reliable Results



## Enhanced Reaction System

- Non-touch vortexer mixing for reaction solution without any contamination
- Precisely controlled temperature for magnetic separation, incubation & reaction by maintenance free direct-solid heating technology
- Reliable clinical performance with enhanced PMT chemiluminescent detection system



## Efficient Washing System

- Interior & exterior probe washing with high pressure pre-heated water
- Programmed enhanced washing with detergent for reagent and sample probes
- 4-phase washing for magnetic separation unit with high efficiency

## Stable Reagent Cartridges Utility

- 2°C~8°C constant cooling compartment for reagents cartridges
- Membrane sealed liquid reagent ready-to-use
- Intelligent magnetic particles mixing by reagent carousel rotation



## Intelligent Indication for Daily Operation

- Real-time indication of cuvettes, waste and wash buffer status
- Real-time monitoring for substrates consumption
- Easy operation for replacing cuvettes, wash buffer and placing new substrates



# Advanced Software

## User-friendly Interface

- Unified platform of BS-2000M series, BS-800M series, BS-480 and future system
- Real-time status monitoring for analytical unit, SDM and carousels
- Intelligent distribution for STAT samples



## Real-time QC Status Monitoring

- Westgard Rules and Two-Control Evaluation check
- Levey-Jennings chart and Twin-Plot chart for review
- Real-time alarm and locating when QC result(s) is out of range
- Auto QC setup capability



## Traceable Test Results

- Reagent, calibrator and control information can be recalled from archive history
- Intuitive software design, easy to trace results



## Reflex and Re-run Function

- Pre-defined reflex and Re-run assays will be performed automatically when preset criteria is met



## Test Summary Function

- Test summary during certain period, including calibration, QC, Sample, valid tests and rerun tests
- Facilitate to computation of total test costs within a defined period
- The summary can be archived into excel files or printed to review and backup



## Step-by-step Maintenance Guide

- Scheduled maintenance and maintenance guide for analytical unit
- Ensure performance reliability and reduce unnecessary service calls
- Error report transferrable to service engineers for immediate troubleshooting which can minimize instrument downtime



## Technical Specifications

### General Information

|                       |  |
|-----------------------|--|
| Throughput            | Up to 240 tests/hour   |
| Measurement Principle | Micron superparamagnetic particles platform with alkaline phosphatase (ALP) labeled reagents and AMPPD substrate |

### Sample Assembly

|                 |   |
|-----------------|---|
| Sample Feeder   | Rack system with Sample Delivery Module<br>Up to 300 onboard samples with 30 racks                |
| Sample Probe    | Steel probe, liquid level detection, clot detection, horizontal and vertical collision protection |
| Sample Volume   | 10 $\mu$ L~200 $\mu$ L, with increment of 1 $\mu$ L   |
| Sample Dilution | Automatic dilution with 1:2 ~ 1:40  |

### Reagent Assembly

|                 |   |
|-----------------|---|
| Reagent Disk    | 36 reagent positions, continuous loading, real-time mixing for magnetic particles                   |
| Cooling System  | 2°C~8°C   |
| Reagent Barcode | Automatic scanning  |
| Reagent Probe   | Steel probe, liquid level detection, bubble detection, horizontal and vertical collision protection |
| Reagent Volume  | 20 $\mu$ L~135 $\mu$ L, with increment of 1 $\mu$ L   |

### Substrate System

|                   |  |
|-------------------|--|
| Substrate Bottle  | 300/500 tests per bottle, bottom piercing    |
| Substrate Loading | Capable of continuous loading during testing |
| Substrate Heating | Constant heating                             |
| Aspiration Volume | 200 $\mu$ L                                  |

### Magnetic Separation Unit

|                     |  |
|---------------------|--|
| Separation Assembly | 4-phase magnetic separation                        |
| Temperature         | Accuracy: 37 $\pm$ 0.3°C, fluctuation: $\pm$ 0.1°C |

### Measurement and Reaction System

|                 |  |
|-----------------|--|
| Detection Mode  | Photon counting                                    |
| PMT Calibration | LED reference module                               |
| Signal Detector | Photomultiplier (PMT) photon counter               |
| Reaction Disk   | 150 positions                                      |
| Temperature     | Accuracy: 37 $\pm$ 0.3°C, Fluctuation: $\pm$ 0.1°C |
| Mixing Unit     | Non-touch vortexer mixing                          |

### Cuvette Loading and Waste Collection

|                  |  |
|------------------|--|
| Cuvette          | Disposable plastic cuvettes  |
| Cuvette Loading  | Trays continuous loading and offloading, 126 cuvettes per tray, up to 9 trays (1134 cuvettes) can be loaded in one batch |
| Waste Collection | 2 waste boxes, overflow detection  |

### Working Conditions

|                       |   |
|-----------------------|---|
| Power Supply          | 110V/115V~, 60Hz; 220V-240V~, 50Hz;<br>220V/230V~, 60Hz   |
| Input Power           | 3200VA  |
| Operating Temperature | 15°C~30°C   |
| Relative Humidity     | 35%RH~80%RH, without condensation   |
| Dimension             | Main Unit: 1450mm*1020mm*1200mm (W*D*H)<br>SDM: 710mm*1020mm*1200mm (W*D*H)<br>Total: 2150mm*1020mm*1200mm (W*D*H)<br>Weight: 750kg |
| Altitude Height       | -400m~2000m<br>(optional external compressor for 2000m~4000m)   |