GeneXpert[®] Infinity-80

5 (E. S.

111111

STEP IN FIFTH

The World's Highest Throughput Molecular Diagnostic System



....

111

11.

1.1

101

i 🖬 🕬



INNOVATION Engineered to fit you

The GeneXpert[®] Infinity-80 is a completely new system. Re-engineered from the ground up for increased capacity, improved efficiency of operation, and even greater reliability - all in the same footprint as the Infinity-48. And the system is designed to work the way you work - batch testing, on-demand or STAT. The Infinity-80 is ready 24/7.

Ergonomic, adjustable touch screen

Fast, efficient sample scanner

Conveyer loading zone handles up to 4 cartridges simultaneously

Module doors eliminate pneumatics for lower maintenance and footprint

Smart waste disposal doesn't interrupt operation

Easily prioritize a STAT sample against the testing queue for the fastest results

許許對對

.

000



GAME CHANGER : CEPHEID GENEXPERT INFINITY-80

Ten years of focused R&D investment is behind the GeneXpert® Infinity-80. With its exclusive 'load and go' technology, the Infinity couples total walkaway operation with complete random access availability. It's the ultimate automated molecular solution for your core laboratory that will have lasting impact on lowering costs and improving patient care throughout your hospital.



The GeneXpert[®] Infinity-80 doesn't need you to keep it company. By automating every aspect, it maximizes the efficiency and resource optimization in your lab. Just scan a disease-specific GeneXpert[®] cartridge into the system, load it, and our Xpertise[™] software takes care of the rest. Runs can be prioritized and time to results identified, all from a single touchscreen. The Infinity is also amazingly flexible. Run single-patient samples as they arrive or in batches. Your choose what's optimal for your workflow, or workload - with the capability of over 2,000 tests per day. The Xpertise software interfaces with any Laboratory Information or Hospital Information System network, delivering real-time answers.



Responsible diagnostic testing demands proven technologies. With more than 10 million Anthrax tests run without a single false positive, no other molecular system has a proven track record that can compete with the GeneXpert System.

FAC1

02:17 20:00 12:40 01:014:00

The United States Postal Service Biohazard Detection System has run over ten million GeneXpert te without a single false positive.



GAME CHANGER : CEPHEID GENEXPERT INFINITY-80





What happens when LEAN principles meet cutting-edge innovation and design? You maximize results per square foot and gain efficiencies that extend beyond the

The GeneXpert Infinity-80 is an investment that pays you back. The modular design and high throughput capability keeps testing costs low while reducing maintenance and maximizing uptime. And the diagnostic accuracy of the Infinity-80 means you'll shrink risk while strengthening the most important bottom line of all: the health of your patients.



For our full GeneXpert family, please visit www.cepheidinternational.com

GAME CHANGER : CEPHEID GENEXPERT INFINITY-80



R

Specifications

Infinity-80

Automated, multi-mode molecular diagnostic analyzer for on-demand access processing of Rt-PCR based clinical diagnostics.

Number of Modules 16, 24, 32, 48, 56, 64, 72, 80

Reagents

Cepheid test-specific, ready-to-use, barcoded kit

Test Throughput (theoretical max)

2,300 in 24 hours

Walk-Away Capacity

183 samples per Infinity (80 loaded in modules, 80 queued, 19 conveyor, 4 shuttle)

QC Method

GeneXpert[®] test includes internal controls to monitor sample, processing and performance.

Sample Barcode Type Code 128, Code 39

Power Requirements

Line 200 - 240 VAC, 20A Voltage: Frequency: 50 Hz – 60 Hz +/- 0.5%

Power Consumption

Max Current 16 Amps @ 200 VAC

10% – 95% RH. non-condensing

15 – 30°C

Heat Output:

Noise Output:

Maximum 530 W/hr or 450 Kcal/hr

< 69 dB

www.Cepheidinternational com

108" (267.2 cm) L 35" (87.9cm) W 78" (200.9 cm) H

2,100 pounds = 954.5 kg

904 Caribbean Drive Sunnyvale, CA 94089 USA +1.888.336.2743

+1.408.541.4191 +1.408.541.4192

Vira Solelh 81470 Maurens-Scopont, France +33.563.82.53.00 FAX +33.563.82.53.01 EMAIL cepheid@cepheideurope.fr

The purchase of this product allows the purchaser to use it for the performance of diagnostic services for human in vitro diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby. No other rights are conveyed expressly, by implication or estoppel to any other patents. Futhermore, no rights for resale are conferred with the purchase of this product.