

STANDARD F 200

- » Convenient and powerful immunoassay analyzer.
- » F200 is a user friendly designed FIA analyzer. Its compact design and convenience features will make your lab-work easier and smoother.



» Technical Specification



Model	STANDARD™ F200
Test method	Fluorescent immunoassay (FIA)
Analysis	Quantitative / Qualitative tests
Test capacity	1 test
Test mode	STANDARD TEST, READ ONLY
Power	AC/DC Adapter
Display	7" Color touch screen
Printer	Built-in
Connectivity	HL7 v2.6(PCD-01) / POCT1-A
Auto-ID	2D Barcode
Accessories	Keyboard / Barcode scanner
Dimension	200 x 240 x 205 mm
Weight	2.5 kg