

STANDARD F 2400

» The best way to reduce turn-around time and improve service quality of your laboratory.



» Technical Specification



Model	STANDARD™ F2400
Test method	Fluorescent immunoassay (FIA)
Analysis	Quantitative / Qualitative tests
Test capacity	70 tests per hour
Test mode	STANDARD TEST
Power	AC/DC Adapter
Display	10" Color touch screen
Printer	Built-in
Connectivity	HL7 v2.6(PCD-01) / POCT1-A
Auto-ID	2D Barcode
Accessories	Keyboard / Barcode scanner
Dimension	510 x 566 x 297 mm
Weight	20.0 kg